

Marketing Authorisation Holder and Manufacturer

Marketing Authorisation Holder:
Mylan Products Ltd. 20 Station Close, Potters Bar, Herts, EN6 1TL, United Kingdom

Manufacturer:
Aspen Notre Dame de Bondeville, 1 rue de l'Abbaye, F-76960 Notre Dame de Bondeville, France.

This leaflet was last revised in November 2023



ART-5428-02

STEP BY STEP GUIDE TO USING ARIXTRA

Types of safety syringe

There are two types of safety syringes used for Arixtra, designed to protect you from needle stick injuries following injection. One type of syringe has an **automatic** needle protection system and the other type has a **manual** needle protection system.

Parts of the syringes:


- ① Needle shield
- ② Plunger
- ③ Finger-grip
- ④ Security sleeve

Picture 1. Syringe with an **automatic** needle protection system




Syringe with a **manual** needle protection system

Picture 2. Syringe with a **manual** needle protection system



Picture 3. Syringe with a **manual** needle protection system showing security sleeve being pulled over needle **AFTER USE**




Instructions for use

These instructions are for both types of syringes (automatic and manual needle protection system). Where the instruction for a syringe is different this is clearly stated.

- 1. Wash your hands thoroughly** with soap and water and dry them with a towel.
- 2. Remove the syringe from the carton and check that:**
 - the expiry date has not passed
 - the solution is clear and colourless to slightly yellow and doesn't contain particles
 - the syringe has not been opened or damaged

- 3. Sit or lie down in a comfortable position.** Choose a place in the lower abdominal (tummy) area, at least 5 cm below your belly button (picture A). **Alternate the left and right side** of the lower abdominal area at each injection. This will help to reduce the discomfort at the injection site.




Picture A

If injecting in the lower abdominal area is not possible, ask your nurse or doctor for advice.

- 4. Clean the injection area with an alcohol wipe.**


- 5. Remove the needle shield,** by first twisting it (picture B1), and then pulling it in a straight line away from the body of the syringe (picture B2). **Discard the needle shield.**



Picture B1


Important note

- **Do not touch the needle** or allow it to touch any surface before the injection.
- It is normal to see a small air bubble in this syringe. **Do not try to remove this air bubble before making the injection** - you may lose some of the medicine if you do.




Picture B2

- 6. Gently pinch the skin that has been cleaned to make a fold.** Hold the fold between the thumb and the forefinger during the entire injection (picture C).




Picture C

- 7. Hold the syringe firmly by the finger grip.** Insert the full length of the needle at right angles into the skin fold (picture D).



Picture D


- 8. Inject ALL of the contents of the syringe by pressing down on the plunger as far as it goes** (picture E).



Picture E

Syringe automatic system

- 9. Release the plunger** and the needle will automatically withdraw from the skin and go back into the security sleeve where it will be locked permanently (picture F).



Picture F

Syringe manual system

- 9. After the injection hold the syringe in one hand by gripping the security sleeve,** use the other hand to hold the finger grip and pull firmly back. This unlocks the sleeve. Slide the sleeve up the body of the syringe until it locks into position over the needle. This is shown in Picture 3 at the beginning of these instructions

Do not dispose of the used syringe in the household waste. Dispose of it as your doctor or pharmacist has instructed.

VIATRIS				Date: 23 NOV 2023	Time: 17:00
Description	arixtra, 10 mg/0.8 ml, 11 mg/0.8 ml, 6 mg/0.4 ml, 7.5 mg/0.6 ml, 5 mg/0.4 ml, 10			No. of colours	3
Component Type	Leaflet/IFU - Folded	Site Barcode/DataMatrix	TBC	Page Count	2 of 2
Affiliate Item Code	ART-5428-02	Viatrix SAP No.	400557889, 400557887, 400557890	Colours	Black Magenta Cyan
Superseded Affiliate Item Code	322802	Vendor Job No.	N.A.	Non-Print Colours	Keyline
TrackWise/GLAMS Job No.	ART-5428-02	Artwork Proof No.	1	Equate CMYK with	
MA No.	PL68 46302/0222, PL68 46302/0223, PL68 46302/0229	Client Market	United Kingdom	Main Font	Frutiger Next
Supplier SAP No.	N.A.	Barcode Info	N.A.	Body Text Size	8 pt
New Supplier Code	TBC	3D Render ID	N.A.	Dimensions	480x200 mm
Superseded Supplier Code	TBC			Min Text Size used	8 pt
Packing Site/Printer	Aspen NDB SAS (Notre-Dame-de-Bondeville - FR)			Keyline/Drawing No.	NB-NOT-011-E-FM_ASPEN_01_LAYOUT1
Sign-offs					